

PUBLIC EXHIBIT 20

Shop

Health Features

Experience

For Organizations

There are no items in your cart

THE PULSE *Blog*



Author: [Oura Team](#) | [October 3, 2024](#)

Oura Ring 4 is our most personal, accurate, and comfortable smart ring yet.

With scientifically validated sensing capabilities, Oura Ring 4 is built around [Smart Sensing](#), a breakthrough technology that adapts to your unique physiology to deliver accurate, continuous data.

Featuring recessed sensors and a fully titanium inner and outer coating, Oura Ring 4 is sleeker and more comfortable than ever to wear all day long, from workout to night out. In keeping with our mission to

BENEFITS BEFORE THE END OF THE YEAR!


HIGH BLOOD PRESSURE IN PREGNANCY: WHAT YOU NEED TO KNOW

6 LIFESTYLE HABITS TO MANAGE BLOOD PRESSURE NATURALLY

WHAT IS HYPERTENSION? SPOT THE SIGNS, LOWER YOUR RISK & PROTECT YOUR HEART

Join the Community

Shop Now



make health a daily practice for every body, Oura Ring 4 is also available in a wider range of sizes—now extending from size 4 to size 15.

~~Sleeker. Smarter. Made for you.~~

[Shop Oura Ring 4](#)

Furthermore, all Oura members with an Oura Ring Gen3 or Oura Ring 4 will experience upgrades to several features, including Daytime Stress, Automatic Activity Detection, and Cycle Insights, as well as a newly [redesigned Oura App](#).

“Oura Ring is designed to be a personal health companion, giving our members a holistic view of their health by pairing personal data with meaningful, actionable insights to facilitate long-term changes,” says Tom Hale, Oura CEO. “The next-generation Oura Ring 4, reimagined Oura App, and the development of the Smart Sensing platform represent significant steps forward for wearables across accuracy, comfort, personalization, and design.”

Oura Ring 4 starts at \$349, with Oura Membership pricing at \$5.99 per month or \$69.99 per year.

Keep reading to learn more.

Oura Ring 4: More Accurate Than Ever Before



Crafted with precision, the [accuracy of Oura Ring](#) has been validated by research institutions across the globe, such as the University of Tokyo and the [National University of Singapore](#). In Oura Ring 4, the new Smart Sensing technology provides an unparalleled degree of accuracy and personalization.

This breakthrough health sensing platform works with the ring's research-grade sensors to learn your unique physiology. The multiple sensors inside the ring form an 18-path multi-wavelength photoplethysmography (PPG) subsystem to provide highly accurate, continuous biometric data, including resting heart rate, [sleep stages](#), [heart rate variability \(HRV\)](#), temperature trends, and more, throughout the day and night.

LEARN MORE: [Smart Sensing: Unparalleled Accuracy and Personalization in Oura Ring 4](#)

A Wearable That's Truly Wearable, Day and Night



Oura Ring 4 is sleeker and more comfortable than its predecessors, allowing it to blend seamlessly into your lifestyle. With a fully titanium coating inside and out, it's gentle to the touch, even on sensitive skin, and slips on and off your finger with ease.

Plus, the powerful health sensors are now recessed, reducing friction and pressure on your fingers when shaking hands or gripping objects. With the lower profile and advancements in comfort of Oura Ring 4, you might even forget you're wearing it—until the compliments start coming.

Choose from Expanded Sizes & Stunning Finishes



Find the right fit for you *and* your personal style. Oura Ring 4 comes in an expanded range of sizes (4-15), the widest size range available on the smart ring market. Select from six gorgeous finishes, including a new, deep black that's more durable than before.

LEARN MORE: [Inside the Ring: The Making of Oura Ring 4](#)

Enjoy Up to 8 Days of Battery Life



Never miss a beat. Oura Ring 4 delivers up to eight days of battery life on just a single charge, and fully charges in up to 120 minutes, depending on current battery level.*

Plus, the charger itself has been redesigned, with an aluminum body and grippy mat on the bottom that provides more stability and improved anchoring as it charges your ring.


Experience a Redesigned Oura App



Currently rolling out to all members, the award-winning Oura App has been redesigned to bring your short- and long-term health into focus, with three new, streamlined tabs: Today, Vitals, and My Health. This easier-to-navigate setup provides an at-a-glance view into your holistic wellness data on a daily basis and trends over time.

READ MORE: [Introducing the New Oura App Design](#)

New Feature Upgrades



Over time, we've significantly invested in improving our software experience. As of October 3, 2024, new software updates will begin rolling out to both Oura Ring 4 and Gen3 members.

[Automatic Activity Detection \(AAD\)](#) now includes **average heart rate and heart rate zones for 40+ activities**, automatically giving you a more comprehensive view of your movement throughout the day.

Cycle Insights now includes **Fertile Window**, a new feature that gives women who are trying to conceive even more information about their cycle and their fertility journey.

Daytime Stress now provides more context about the relationship between your daytime stress levels and your tags, movement, and activities. This will help you understand how your daily habits affect your stress and your recovery.

Oura Ring 4 vs. Oura Ring Gen3

How does Oura Ring 4 compare to Oura Ring Gen3? The comparison chart below provides specifics, such as the difference in material and design, the sensing platform, and sizing.

In addition, Oura Ring 4 is now capable of storing up to two weeks of data, depending on the type and frequency of use.

Compare Oura Ring Models

Style	 Oura Ring 4	 Oura Ring Gen3
Design & Material	Fully titanium exterior & interior Recessed interior sensors at 0.3mm	Titanium exterior, epoxy interior Raised interior domes (sensors) at 1.3mm
Sensing Platform	Accuracy is powered by Smart Sensing that dynamically selects the optimal signal path	Measurements powered by a fixed, static configuration of optical paths
Battery	Up to 8 days	Up to 7 days
Sizing	Sizes 4-15	Sizes 6-13
Membership Pricing	\$5.99/month	\$5.99/month

**All battery claims depend on Oura Ring size, Oura App settings, feature configuration, activity, age of battery, usage, and many other factors; actual results will vary. Battery life and charge cycles vary by use and settings. Battery tests were conducted across*

various sizes and scenarios. See [Ring Battery tips](#) for more information.

Featured Articles



COMMUNITY
MAXIMIZE YOUR
FSA/HSA
BENEFITS
BEFORE
THE END
OF THE
YEAR!



HEALTH
OURA
CYCLE
INSIGHTS
UPDATE:
FASTER,
MORE
ACCURATE
PERIOD
&
OVULATION
TRACKING



HEALTH
CYCLE
INSIGHTS: KNOW
YOUR
CYCLE
PHASES,
FERTILE
WINDOW
, PERIOD
PREDICTION, AND
MORE



NEWS
OURA
FILES
ITC
ACTION
AGAINST
SAMSUNG,
REEBOK
(REEBOK
SMART
RING),
ZEPP
HEALTH
(AMAZFIT RING),
AND
NEXXBAND
FOR
PATENT
INFRINGEMENT



Instant Checkout



Major Cards Accepted

Receive articles, tips, and offers from Oura

Email Address

We care about protecting your data.

Read more in our [Privacy Policy](#).

We care about protecting your data. Read more in our [Privacy Policy](#).

Our Company

About Us

Leadership

Medical Advisory Board

Careers

Newsroom

Support

Help Center

Sizing

Recycling Program

Flexible Spending

Heart Rate Monitoring

My Account

Oura on the Web

Contact

Partner With Us

For Organizations

Partnerships

Developers

Connect

The Pulse *Blog*

Facebook

Instagram

Pinterest

TikTok

X

YouTube

[Terms & Conditions](#)

[Privacy Policy](#)

[Accessibility](#)


[FCC Compliance](#)

[IP Notice](#)

[Security Center](#)

© 2025 Ōura Health Oy. All Rights Reserved.

ŌURA, OURA, OURA RING, and Ō are trademarks of Oura Health Oy and may not be used without permission.



Shop

Health Features

Experience

For Organizations

There are no items in your cart

THE PULSE *Blog*



Author: Emilina Lomas | March 18, 2025

Expert review by: Xi Zhang, PhD, VP, Health Sensing

[Oura Ring](#) has become a leader in personal health tracking with its advanced, discreet technology. It goes beyond basic step counting, delving deep into [sleep](#), [cardiovascular health](#), [cycle tracking](#), [stress](#), [illness detection](#), and more than 30 other biometrics.

[Oura Ring 4](#) represents the pinnacle of wearable health technology, integrating and minimizing advanced sensors, proprietary algorithms, and continuous data analysis into a comfortable ring form-factor to deliver highly accurate health insights.

BENEFITS BEFORE THE END OF THE YEAR!


HIGH BLOOD PRESSURE IN PREGNANCY: WHAT YOU NEED TO KNOW

6 LIFESTYLE HABITS TO MANAGE BLOOD PRESSURE NATURALLY

WHAT IS HYPERTENSION? SPOT THE SIGNS, LOWER YOUR RISK & PROTECT YOUR HEART

Join the
Community

Shop
Now



Learn more about the cutting-edge technology that powers Oura Ring 4—which, in turn, empowers Oura members to take control of their health and wellbeing.

~~Sleeker. Smarter. Made for you.~~

[Shop Oura Ring 4](#)

Advanced Sensor Suite

Oura Ring 4 features a completely redesigned multi-sensor system, offering unparalleled accuracy.

- **18-path multi-wavelength photoplethysmography (PPG):**
 - Red and infrared LEDs measure blood oxygen levels (SpO2) while you sleep.
 - Green and infrared LEDs track heart rate (HR) and heart rate variability (HRV) 24/7, as well as respiration rate during sleep.
- **Digital temperature sensor:** measures skin temperature variations for [early signs of sickness](#), [menstrual cycle tracking](#), and [recovery insights](#).
- **3D accelerometer:** monitors [movement and](#)

[activity 24/7](#), capturing steps, exercise intensity, and nighttime restlessness.

- **Groundbreaking [Smart Sensing](#) technology:** dynamically adjusts PPG sensor configurations based on multiple contexts and conditions, significantly improving signal fidelity and battery efficiency.



Proprietary Algorithms for Data Analysis

The research-grade sensors power state-of-the-art machine-learning algorithms to deliver the most accurate, continuous health metrics and insights, day and night—all while optimizing battery life.

- Algorithms adapt to your personal baseline, adjusting contributors to your [Sleep](#), [Activity](#), and [Readiness Scores](#) over time.
- Oura's advanced [machine-learning-based sleep staging algorithm](#) has been [scientifically validated](#) against polysomnography (PSG)—the gold

standard of sleep tracking—allowing it to generate precise sleep stage graphs (hypnograms).

- [Automatic Activity Detection \(AAD\)](#) intelligently detects 40+ workouts without manual input, delivering average heart rate and heart rate zones after the activity.
- Smart Sensing reduces data gaps, leading to a 120% improvement in signal quality for [blood oxygen sensing \(SpO2\)](#), 31% improvement in [nighttime heart rate tracking](#), and a 7% improvement in daytime heart rate accuracy.

RELATED: [Technology in Oura Ring Gen3](#)

Bluetooth Connectivity

The bluetooth technology in Oura Ring is only active during a small portion (well below 1%) of each day. Data is transmitted more continuously when the ring syncs with the app and during firmware updates. [Airplane mode](#) disables all radio transmission from Oura Ring.

The limit for SAR (Specific Absorption Rate) for a device such as Oura Ring is 2.0 W/kg for head and body, and Oura's SAR level is 0.0003 W/kg. In other words, the electromagnetic field (EMF) levels of Oura Ring are about *3,000 times* lower than a typical smartphone.

Dimensions and Weight

- Width: 7.90 mm
- Thickness: 2.88mm
-

Weight Range: 3.3 — 5.2 g for full range of sizes 4 — 15

Water Resistance

- Water resistant to 100m/328ft

Battery


- Up to 8 days of battery life*
- Typically charges fully in 20-80 minutes, depending on the battery level
- Lithium-ion rechargeable battery

**All battery claims depend on Oura Ring size, Oura App settings, feature configuration, activity, age of battery, usage, and many other factors; actual results will vary. Battery life and charge cycles vary by use and settings. Battery tests were conducted across various sizes and scenarios.*


Featured Articles



HEALTH
OURA
CYCLE
INSIGHT
S
UPDATE:
FASTER,



HEALTH
CYCLE
INSIGHT
S: KNOW
YOUR
CYCLE
PHASES,



MEET OURA
INSIDE
THE
RING:
WHAT
FUELS
YOUR



MEET OURA
INTRODU
CING
CUMULA
TIVE
STRESS:
UNDERS

MORE
ACCURATE
PERIOD
&
OVULATION
TRACKING

FERTILE
WINDOW
, PERIOD
PREDICTION,
AND
MORE

NUTRITION
LEVEL

TANDING
THE
IMPACT
OF
LONG-
TERM
STRESS

Instant Checkout

Major Cards Accepted

Receive articles, tips, and offers from Oura

Email Address

We care about protecting your data.
Read more in our [Privacy Policy](#).
We care about protecting your data. Read more in our [Privacy Policy](#).

Our Company

- About Us
- Leadership
- Medical Advisory Board
- Careers
- Newsroom

Support

- Help Center
- Sizing
- Recycling Program
- Flexible Spending
- Heart Rate Monitoring
- My Account
- Oura on the Web
- Contact

Partner With Us

For Organizations

Partnerships

Developers

Connect

The Pulse *Blog*

Facebook

Instagram

Pinterest

TikTok

X

YouTube

Terms & Conditions

Privacy Policy

Accessibility


FCC Compliance


IP Notice

Security Center

© 2025 Ōura Health Oy. All Rights Reserved.

ŌURA, OURA, OURA RING, and Ō are trademarks of Oura Health Oy and may not be used without permission.





Shop

[Oura Help](#) > [Oura Ring 4](#) > [Rings](#)

Articles in this section



Oura Ring 4

Oura Ring 4 is the next evolution of the revolutionary smart ring, designed to help you get up close and personal with your health. Oura Ring 4 is built around a breakthrough sensing technology, Smart Sensing, which adapts to your finger to deliver accurate, continuous data, day and night.

- Oura Ring 4 is available through the [Oura online store](#) and various [authorized retailers](#)
- Oura Ring 4 Ceramic:
 - All colors (Midnight, Cloud, Petal, and Tide) are available through the [Oura online store](#) and [Best Buy US](#) and [Best Buy Canada](#)
 - Cloud, Petal, and Tide are available at various [authorized retailers](#)

If you're setting up your Oura Ring 4 as a new or current member, see the instructions in [Set Up an Oura Ring](#).

Also see [Oura Ring 4 FAQs](#).

[Hardware](#)

[Software](#)

[Oura Ring 4 vs. Oura Ring 4 Ceramic](#)

[More Information](#)

Hardware

Oura Ring 4 features a fully titanium design and recessed sensors on the interior, enhancing both durability and wearability. Oura Ring 4 is available in six finishes: Black, Silver, Stealth, Brushed Silver, Gold, and Rose Gold.

Oura Ring 4 Ceramic is made from high-performance zirconia ceramic, a premium material selected for its scratch resistance, luster, and color. Oura Ring 4 Ceramic is available in four colors: Midnight, Cloud, Petal, and Tide.

All Oura Ring 4 models are available in sizes 4-15.

The sizing for Oura Ring 4 is slightly different from Oura Ring Gen3. Make sure to use the right sizing kit for the most accurate fit.

- Oura Ring 4 Sizing Kit: white ring sizers, sizes 4–15
- Oura Ring Gen3 Sizing Kit: black ring sizers, sizes 6–13

For the best fit and most accurate readings, align your ring so that the sensors are against the palm side of your finger.

Sensors

- Smart Sensing is a sensing platform designed to give you the most complete picture of your health. The research-grade sensors work with an algorithm that adapts to the unique structure and skin tone of your finger and your current activity to deliver the most accurate, continuous data, day and night—all while optimizing battery life
- Red and infrared LEDs measure blood oxygen while you sleep
- Green and infrared PPG sensors alternate to measure heart rate 24/7 and heart rate variability and respiration rate during your sleep
- Digital temperature sensor measures average body temperature variations
- Accelerometer tracks movement and activity 24/7

External Surface Materials

- Oura Ring 4:
 - Titanium with physical vapor deposition (PVD) coating in Black, Silver, Brushed

Silver, Gold, and Rose Gold

- Titanium with Diamond-like Carbon (DLC) coating in Stealth
- Oura Ring 4 Ceramic:
 - Zirconia ceramic

Inner Surface Materials

- Non-allergenic, titanium inner shell. This is a medical grade, BPA-free inner surface
- Shock resistance and bonding between printed circuit board and ring outer frame

Connectivity

- Bluetooth Low Energy (Bluetooth Smart®)
- EMF-Safe and allows Airplane Mode
- FCC certified

Other Specifications

- Battery: Typically 5-8 days. Battery life and longevity depend on Oura Ring size, Oura App settings, feature configuration, activity, battery age, usage, and many other factors; actual results will vary. Battery life and charge cycles vary by use and settings. Battery tests were conducted across various sizes and scenarios. See Oura Ring 4 Battery Life and Ring Battery Tips for more information
- Charging time: Typically between 20-80 minutes depending on the ring starting charge level
- Water Resistance: 100m/328ft, suitable for sauna and water sports but not diving
- Operating Temperature: -10–54 °C / 14–129 °F. Charging at room temperature recommended

See more information in our [Product Safety & Use article](#).

Software

Oura Ring 4 with active membership includes:

- [Detailed sleep analysis](#)

- 24/7 heart rate tracking (Daytime, Nighttime, Activity)
- Integration with third-party health and wellness apps like Strava
- Advanced temperature monitoring
- Blood Oxygen Sensing (SpO2)
- Daytime Stress
- Dynamic activity goals with Automatic Activity Detection
- Heart Health (Cardiovascular Age, Cardio Capacity (VO2 Max))
- Metabolic Health (Glucose Tracking, Meals)
- Oura Labs
- Resilience
- Weekly, monthly, quarterly, yearly, and anniversary reports
- Women's Health (Cycle Insights, Pregnancy Insights)

Body clock and chronotype, Explore content, and personalized activity goals are also available for Gen2 and Gen3 rings.

Regardless of the type of membership you have, your Oura Membership will continue to work seamlessly with Oura Ring 4. Since your membership is tied to your Oura account (not the ring itself), all of your data and benefits will carry over. Once you pair your new Oura Ring 4, you'll see all your data in the Oura App.

See [Oura Membership](#) for more information and pricing.

Oura Ring 4 vs. Oura Ring 4 Ceramic

	Oura Ring 4	Oura Ring 4 Ceramic
Design & Material	• Fully titanium exterior • Titanium interior • Recessed interior sensors at 0.3mm	• Ceramic exterior • Titanium interior • Recessed interior sensors at 0.3mm
	• Black • Brushed Silver	• Cloud

Finishes/Colors	• Gold • Silver • Stealth • Rose Gold	• Midnight • Petal • Tide
Sensing Platform	Accuracy is powered by Smart Sensing that dynamically selects the optimal signal path	Accuracy is powered by Smart Sensing that dynamically selects the optimal signal path
Sensors	• Red and infrared LEDs measure blood oxygen levels • Green and infrared LEDs alternate to measure heart rate and heart rate variability 24/7 and respiration rate • Digital sensor measures temperature trends • Accelerometer tracks movement and activity 24/7	• Red and infrared LEDs measure blood oxygen levels • Green and infrared LEDs alternate to measure heart rate and heart rate variability 24/7 and respiration rate • Digital sensor measures temperature trends • Accelerometer tracks movement and activity 24/7
Dimensions	• Width: 7.9mm • Thickness: ~2.88mm • Weight: 3.3-5.2 grams (depending on ring size)	• Width: 7.9mm • Thickness: 3.51mm • Weight: 5.1 - 8.1 grams (depending on ring size)

More Information

Product Safety & Use

Oura Membership

Previous Oura Ring Generations

Oura Ring Finishes

Ring Battery Tips

Was this article helpful?

☐ Yes

☐ No

1134 out of 1719 found this helpful

Return to top ☐

Recently viewed articles

Blood Oxygen Sensing (SpO2)

Nap Detection

Sleep Stages

Live Heart Rate

Product Safety & Use

Related articles

Oura Membership

How to Choose the Right Oura Ring Size

Retail Locations

Ring Battery Tips

Oura Ring Generation 3